

# CATALOG DIGITAL IMAGING SUPPLIES

- BARRIER ENVELOPE
- CARTON PROTECTIVE COVERS
- PHOSPHOR PLATES

## Barrier Envelopes

- Premium Quality, Soft to touch.
- Reduce Cross-Infection.
- Keep the Phosphor Storage Plate Clean.
- Easy to use.
- Low cost, Disposable Use.
- Fit all brands!



**500 pcs / box**  
 Available Sizes:  
 Size #2, 41 x 31 mm  
 Size #1, 40 x 24 mm  
 Size #0, 35 x 22 mm



## Carton Protective Covers

- Protects the Phosphor Storage Plates from mechanical damage and information loss caused by ambient light.
- Easy to use.
- Reduce Cross-Contamination.
- Fit most brands - Scan-X, Apixia, Vista Scan, Gendex, Digora and more!




Available Sizes:  
 Size #2, 41 x 31 mm



## Phosphor Plates

- Reusable - The Phosphor Plates are re-usable 2000 times to maximize your investment.
- High Resolution - With our Phosphor Plates you will see a much clearer image!
- Affordable - Premium quality for any budget!



Available Sizes:  
 Size #2, 41 x 31 mm  
 Size #1, 40 x 24 mm  
 Size #0, 35 x 22 mm  
 Size #3, 54 x 27 mm



